

For Research Use Only

AFF4 Polyclonal antibody

Catalog Number: 13197-1-AP



Basic Information

Catalog Number:

13197-1-AP

Size:

150ul , Concentration: 650 ug/ml by Nanodrop and 347 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3931

GenBank Accession Number:

BC025700

GeneID (NCBI):

27125

UNIPROT ID:

Q9UHB7

Full Name:

AF4/FMR2 family, member 4

Calculated MW:

365 aa, 40 kDa

Observed MW:

150 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : Jurkat cells, HL-60 cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Jurkat cells were subjected to SDS PAGE followed by western blot with 13197-1-AP (AFF4 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.