

For Research Use Only

PI3K p101 Polyclonal antibody

Catalog Number: 13329-1-AP **9 Publications**



Basic Information

Catalog Number:

13329-1-AP

Size:

150ul, Concentration: 200 µg/ml by Nanodrop and 147 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4137

GenBank Accession Number:

BC028212

GeneID (NCBI):

23533

UNIPROT ID:

Q8WYR1

Full Name:

phosphoinositide-3-kinase, regulatory subunit 5

Calculated MW:

880 aa, 97 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse, rabbit

Notable Publications

Author	Pubmed ID	Journal	Application
Yiyi Jin	29285087	Exp Ther Med	WB
Xiaoqin Zhu	34887935	Evid Based Complement Alternat Med	WB
Hong-Yan Lin	32175001	Chin Med	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

