

For Research Use Only

PI3K p101 Polyclonal antibody

Catalog Number: 13329-1-AP **8 Publications**



Basic Information

Catalog Number: 13329-1-AP	GenBank Accession Number: BC028212	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 200 µg/ml by Nanodrop and 147 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 23533	
Source: Rabbit	Full Name: phosphoinositide-3-kinase, regulatory subunit 5	
Isotype: IgG	Calculated MW: 880 aa, 97 kDa	
Immunogen Catalog Number: AG4137		

Applications

Tested Applications:
ELISA

Cited Applications:
WB

Species Specificity:
human, mouse, rat

Cited Species:
human, mouse, rabbit

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Yiyi Jin	29285087	Exp Ther Med	WB
Xiaoqin Zhu	34887935	Evid Based Complement Alternat Med	WB
Hong-Yan Lin	32175001	Chin Med	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data