

For Research Use Only

LENG8 Polyclonal antibody

Catalog Number: 13336-1-AP

1 Publications



Basic Information

Catalog Number:

13336-1-AP

Size:

150ul , Concentration: 507 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4153

GenBank Accession Number:

BC028048

GeneID (NCBI):

114823

UNIPROT ID:

Q96PV6

Full Name:

leukocyte receptor cluster (LRC) member 8

Calculated MW:

800 aa, 88 kDa

Observed MW:

86-88 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, FC (Intra), ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

WB : A2780 cells, Jurkat cells

Notable Publications

Author	Pubmed ID	Journal	Application
Xiaojuan Ye	26884854	Int J Clin Exp Pathol	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



A2780 cells were subjected to SDS PAGE followed by western blot with 13336-1-AP (LENG8 antibody) at dilution of 1:1200 incubated at room temperature for 1.5 hours.



Jurkat cells were subjected to SDS PAGE followed by western blot with 13336-1-AP (LENG8 antibody) at dilution of 1:2500 incubated at room temperature for 1.5 hours.



1x10⁶ HEK-293 cells were intracellularly stained with 0.25 ug LENG8 Polyclonal antibody (13336-1-AP) and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2)(red), or 0.25 ug Isotype Control (blue). Cells were fixed and permeabilized with Transcription Factor Staining Buffer Kit (PF00011).