

RPP14 Polyclonal antibody

Catalog Number: 13339-1-AP

Basic Information

Catalog Number:

13339-1-AP

Size:

150ul , Concentration: 187 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4159

GenBank Accession Number:

BC012017

GeneID (NCBI):

11102

UNIPROT ID:

O95059

Full Name:

ribonuclease P/MRP 14kDa subunit

Calculated MW:

124 aa, 14 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

RPP14 is a subunit of human ribonuclease P. RPP14 may be part of an exosome activity. This antibody is specific to RPP14.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

