

For Research Use Only

# KCNRG Polyclonal antibody, PBS Only

Catalog Number: 13344-1-PBS



## Basic Information

**Catalog Number:**

13344-1-PBS

**Size:**

100ug, Concentration: 1 mg/ml by Nanodrop;

**Source:**

Rabbit

**Isotype:**

IgG

**Immunogen Catalog Number:**

AG4176

**GenBank Accession Number:**

BC032343

**GeneID (NCBI):**

283518

**UNIPROT ID:**

Q8N5I3

**Full Name:**

potassium channel regulator

**Calculated MW:**

272 aa, 31 kDa

**Observed MW:**

31 kDa

**Purification Method:**

Antigen affinity purification

## Applications

**Tested Applications:**

WB, Indirect ELISA

**Species Specificity:**

human, mouse, rat

## Storage

**Storage:**

Store at -80°C.

**Storage Buffer:**

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 13344-1-AP (KCNRG antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 13344-1-PBS in a different storage buffer formulation.