

For Research Use Only

JAK3 Polyclonal antibody

Catalog Number: 13551-1-AP

2 Publications



Basic Information

Catalog Number:

13551-1-AP

Size:

150ul , Concentration: 67 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4298

GenBank Accession Number:

BC028068

GeneID (NCBI):

3718

UNIPROT ID:

P52333

Full Name:

Janus kinase 3 (a protein tyrosine kinase, leukocyte)

Calculated MW:

1124 aa, 125 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse

Notable Publications

Author	Pubmed ID	Journal	Application
Qiqiang Zhang	33228635	BMC Complement Med Ther	WB
Xue Bai	38464791	J Pharm Anal	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

