

For Research Use Only

# MAK Polyclonal antibody

Catalog Number: 13638-1-AP

1 Publications



## Basic Information

### Catalog Number:

13638-1-AP

### Size:

150ul , Concentration: 400 ug/ml by Nanodrop and 193 ug/ml by Bradford method using BSA as the standard;

### Source:

Rabbit

### Isotype:

IgG

### Immunogen Catalog Number:

AG4571

### GenBank Accession Number:

BC039825

### GeneID (NCBI):

4117

### UNIPROT ID:

P20794

### Full Name:

male germ cell-associated kinase

### Calculated MW:

623 aa, 71 kDa

### Observed MW:

60-62 kDa

### Purification Method:

Antigen affinity purification

### Recommended Dilutions:

WB 1:500-1:1000

## Applications

### Tested Applications:

WB, ELISA

### Cited Applications:

WB

### Species Specificity:

human, mouse, rat

### Positive Controls:

WB : PC-3 cells, HeLa cells

## Notable Publications

Author	Pubmed ID	Journal	Application
Ran Li	28073914	J Biol Chem	WB

## Storage

### Storage:

Store at -20°C. Stable for one year after shipment.

### Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



PC-3 cells were subjected to SDS PAGE followed by western blot with 13638-1-AP (MAK antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.