

For Research Use Only

Glypican 5 Polyclonal antibody

Catalog Number: 13678-1-AP



Basic Information

Catalog Number: 13678-1-AP	GenBank Accession Number: BC039730	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 133 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 2262	
Source: Rabbit	Full Name: glypican 5	
Isotype: IgG	Calculated MW: 572 aa, 64 kDa	
Immunogen Catalog Number: AG4493	Observed MW: 64 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse, rat

Background Information

GPC5, also named as Glypican-5 and Secreted glypican-5, is cell surface proteoglycan that bears heparan sulfate. It is primarily expressed in the brain. 65kd band is Glycosylation; 45kd band is mature form.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data