

For Research Use Only

# TAF8 Polyclonal antibody

Catalog Number: 13693-1-AP



## Basic Information

**Catalog Number:**  
13693-1-AP

**Size:**  
150ul , Concentration: 180 µg/ml by  
Bradford method using BSA as the  
standard;

**Source:**  
Rabbit

**Isotype:**  
IgG

**Immunogen Catalog Number:**  
AG4577

**GenBank Accession Number:**  
BC033728

**GeneID (NCBI):**  
129685

**Full Name:**  
TAF8 RNA polymerase II, TATA box  
binding protein (TBP)-associated  
factor, 43kDa

**Calculated MW:**  
310 aa, 43 kDa

**Purification Method:**  
Antigen affinity purification

## Applications

**Tested Applications:**  
ELISA

**Species Specificity:**  
human, mouse, rat

## Background Information

### Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data