

For Research Use Only

ARHGAP20 Polyclonal antibody

Catalog Number: 13775-1-AP



Basic Information

Catalog Number: 13775-1-AP	GenBank Accession Number: BC039340	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 187 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 57569	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	Full Name: Rho GTPase activating protein 20	
Isotype: IgG	Calculated MW: 133 kDa	
Immunogen Catalog Number: AG4764		

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : human brain tissue,
Species Specificity: human	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human brain tissue were subjected to SDS PAGE followed by western blot with 13775-1-AP (ARHGAP20 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.