

For Research Use Only

# USP6NL Polyclonal antibody

Catalog Number: 13896-1-AP



---

## Basic Information

|  |  |  |
|--|--|--|
| <b>Catalog Number:</b><br>13896-1-AP   | <b>GenBank Accession Number:</b><br>BC042943 | <b>Purification Method:</b><br>Antigen affinity purification |
| <b>Size:</b><br>150ul , Concentration: 333 µg/ml by Bradford method using BSA as the standard; | <b>GeneID (NCBI):</b><br>9712                |  |
| <b>Source:</b><br>Rabbit   | <b>Full Name:</b><br>USP6 N-terminal like    |  |
| <b>Isotype:</b><br>IgG   | <b>Calculated MW:</b><br>828 aa, 94 kDa      |  |
| <b>Immunogen Catalog Number:</b><br>AG4858   |  |  |

---

## Applications

**Tested Applications:**  
ELISA

**Species Specificity:**  
human

---

## Background Information

---

### Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

**\*\*\* 20ul sizes contain 0.1% BSA**

---

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

---

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

---

## Selected Validation Data