

For Research Use Only

JIP2 Polyclonal antibody

Catalog Number: 13912-1-AP



Basic Information

Catalog Number: 13912-1-AP	GenBank Accession Number: BC047527	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 193 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 23542	
Source: Rabbit	Full Name: mitogen-activated protein kinase 8 interacting protein 2	
Isotype: IgG	Calculated MW: 824 aa, 88 kDa	
Immunogen Catalog Number: AG4902		

Applications

Tested Applications:
IHC, ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.


Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human testis using 13912-1-AP (JIP2 antibody) at dilution of 1:50 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human testis using 13912-1-AP (JIP2 antibody) at dilution of 1:50 (under 40x lens).