

For Research Use Only

FZD6 Polyclonal antibody

Catalog Number: 13982-1-AP



Basic Information

Catalog Number:

13982-1-AP

Size:

150ul , Concentration: 193 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5027

GenBank Accession Number:

BC060836

GeneID (NCBI):

8323

UNIPROT ID:

O60353

Full Name:

frizzled homolog 6 (Drosophila)

Calculated MW:

79 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IHC, ELISA

Species Specificity:

human

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human breast cancer using 13982-1-AP (FZD6 antibody) at dilution of 1:50 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human breast cancer using 13982-1-AP (FZD6 antibody) at dilution of 1:50 (under 40x lens).