

For Research Use Only

TRAK1 Polyclonal antibody, PBS Only

Catalog Number: 13987-1-PBS



Basic Information

Catalog Number:

13987-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5055

GenBank Accession Number:

BC015922

GeneID (NCBI):

22906

UNIPROT ID:

Q9UPV9

Full Name:

trafficking protein, kinesin binding 1

Calculated MW:

953 aa, 106 kDa

Observed MW:

64 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 13987-1-AP (TRAK1 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 13987-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human stomach tissue slide using 13987-1-AP (TRAK1 antibody) at dilution of 1:200 (under 10x lens).. This data was developed using the same antibody clone with 13987-1-PBS in a different storage buffer formulation.