

For Research Use Only

DAZ Polyclonal antibody, PBS Only

Catalog Number: 14065-1-PBS



Basic Information

Catalog Number:

14065-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5116

GenBank Accession Number:

BC047617

GeneID (NCBI):

57135

UNIPROT ID:

Q86SG3

Full Name:

deleted in azoospermia 4

Calculated MW:

65 kDa

Observed MW:

60-65 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

DAZ4, encoded by the DAZ gene, is an RNA-binding protein with a role in spermatogenesis. DAZ4 is a testis-specific protein that expression restricted to premeiotic germ cells, particularly in spermatogonia. Defects in DAZ4 may be a cause of spermatogenic failure Y-linked type 2 (SPGFY2). It is a disorder resulting in the absence (azoospermia) or reduction (oligozoospermia) of sperm in the semen, leading to male infertility.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human testis tissue were subjected to SDS PAGE followed by western blot with 14065-1-AP (DAZ4 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 14065-1-PBS in a different storage buffer formulation.