

For Research Use Only

NFRKB Polyclonal antibody, PBS Only

Catalog Number:14086-1-PBS



Basic Information

Catalog Number: 14086-1-PBS	GenBank Accession Number: BC063280	Purification Method: Antigen affinity purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 4798	
Source: Rabbit	UNIPROT ID: Q6P4R8	
Isotype: IgG	Full Name: nuclear factor related to kappaB binding protein	
Immunogen Catalog Number: AG5216	Calculated MW: 1299 aa, 139 kDa	
	Observed MW: 160 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

NFRKB (Nuclear Factor Related to KappaB Binding Protein) plays a role in the recombination of DNA, it helps in the regulation of transcription, DNA-templated, and transcription by RNA polymerase II, and is involved in the regulation of inflammatory responses. It also regulates protein deubiquitination, which is important for protein stability and function.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 14086-1-AP (NFKB antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 14086-1-PBS in a different storage buffer formulation.