

For Research Use Only

IRGC Polyclonal antibody

Catalog Number: 14090-1-AP

1 Publications



Basic Information

Catalog Number:

14090-1-AP

Size:

150ul, Concentration: 240 ug/ml by Nanodrop and 160 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5279

GenBank Accession Number:

BC066939

GeneID (NCBI):

56269

UNIPROT ID:

Q6NXR0

Full Name:

immunity-related GTPase family, cinema

Calculated MW:

50 kDa

Observed MW:

50 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

mouse

Positive Controls:

WB : mouse testis tissue, rat testis tissue

Notable Publications

| Author | Pubmed ID | Journal | Application |
|-------------|-----------|----------|-------------|
| Yuki Kaneda | 35618043 | Dev Biol | WB |

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 14090-1-AP (IRGC antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.