

For Research Use Only

VSIG1 Polyclonal antibody, PBS Only

Catalog Number:14145-1-PBS



Basic Information

Catalog Number:

14145-1-PBS

Size:

100ug , Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5314

GenBank Accession Number:

BC043216

GeneID (NCBI):

340547

UNIPROT ID:

Q86XK7

Full Name:

V-set and immunoglobulin domain
containing 1

Calculated MW:

387 aa, 42 kDa

Observed MW:

65-70 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human stomach tissue were subjected to SDS PAGE followed by western blot with 14145-1-AP (VSI G1 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 14145-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human stomach tissue slide using 14145-1-AP (VSI G1 antibody) at dilution of 1:400 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 14145-1-PBS in a different storage buffer formulation.