

For Research Use Only

# POTEB Polyclonal antibody

Catalog Number: 14169-1-AP

1 Publications



## Basic Information

Catalog Number:

14169-1-AP

Size:

150ul , Concentration: 107 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5368

GenBank Accession Number:

BC066981

GeneID (NCBI):

339010

UNIPROT ID:

Q6S5H4

Full Name:

POTE ankyrin domain family, member B

Calculated MW:

66 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

IHC, ELISA

Cited Applications:

IHC

Species Specificity:

human

Cited Species:

human

## Notable Publications

Author	Pubmed ID	Journal	Application
Redfield Samantha M SM	23841027	Am J Cancer Res	IHC

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

