

For Research Use Only

# SMC5 Polyclonal antibody, PBS Only

Catalog Number: 14178-1-PBS

Featured Product



## Basic Information

Catalog Number:

14178-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5385

GenBank Accession Number:

BC038225

GeneID (NCBI):

23137

UNIPROT ID:

Q8IY18

Full Name:

structural maintenance of chromosomes 5

Calculated MW:

1101 aa, 129 kDa

Observed MW:

120-130 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 14178-1-AP (SMC5 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 14178-1-PBS in a different storage buffer formulation.