

For Research Use Only

# DCLRE1A Polyclonal antibody

Catalog Number: 14210-1-AP



## Basic Information

### Catalog Number:

14210-1-AP

### Size:

150ul , Concentration: 133 µg/ml by Bradford method using BSA as the standard;

### Source:

Rabbit

### Isotype:

IgG

### Immunogen Catalog Number:

AG5425

### GenBank Accession Number:

BC013124

### GeneID (NCBI):

9937

### UNIPROT ID:

Q6PJP8

### Full Name:

DNA cross-link repair 1A (PSO2 homolog, *S. cerevisiae*)

### Calculated MW:

116 kDa

### Purification Method:

Antigen affinity purification

## Applications

### Tested Applications:

ELISA

### Species Specificity:

human, mouse, rat

## Storage

### Storage:

Store at -20°C. Stable for one year after shipment.

### Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

