

For Research Use Only

PARG Polyclonal antibody

Catalog Number: 14227-1-AP

1 Publications



Basic Information

Catalog Number:

14227-1-AP

Size:

150ul , Concentration: 153 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5466

GenBank Accession Number:

BC052966

GeneID (NCBI):

8505

UNIPROT ID:

Q86W56

Full Name:

poly (ADP-ribose) glycohydrolase

Calculated MW:

111 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Notable Publications

Author	Pubmed ID	Journal	Application
Sho Watanabe	36714847	PNAS Nexus	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

