

For Research Use Only

# NEUROG2 Polyclonal antibody

Catalog Number: 14264-1-AP



## Basic Information

|  |  |  |
|--|--|--|
| <b>Catalog Number:</b><br>14264-1-AP   | <b>GenBank Accession Number:</b><br>BC036847 | <b>Purification Method:</b><br>Antigen affinity purification |
| <b>Size:</b><br>150ul , Concentration: 260 µg/ml by<br>Bradford method using BSA as the<br>standard; | <b>GeneID (NCBI):</b><br>63973               |  |
| <b>Source:</b><br>Rabbit   | <b>Full Name:</b><br>neurogenin 2            |  |
| <b>Isotype:</b><br>IgG   | <b>Calculated MW:</b><br>29 kDa              |  |
| <b>Immunogen Catalog Number:</b><br>AG5552   |  |  |

## Applications

**Tested Applications:**  
ELISA

**Species Specificity:**  
human, mouse, rat

## Background Information

### Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data