

For Research Use Only

DPP6 Polyclonal antibody

Catalog Number: 14285-1-AP



Basic Information

Catalog Number:

14285-1-AP

Size:

150ul , Concentration: 450 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5584

GenBank Accession Number:

BC035912

GeneID (NCBI):

1804

UNIPROT ID:

P42658

Full Name:

dipeptidyl-peptidase 6

Calculated MW:

98 kDa

Observed MW:

98-100 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : rat brain tissue, mouse brain tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 14285-1-AP (DPP6 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.