

For Research Use Only

BIK Polyclonal antibody

Catalog Number: 14389-1-AP



Basic Information

Catalog Number:

14389-1-AP

Size:

150ul

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5859

GenBank Accession Number:

BC001599

GeneID (NCBI):

638

UNIPROT ID:

Q13323

Full Name:

BCL2-interacting killer (apoptosis-inducing)

Calculated MW:

18 kDa

Purification Method:

Antigen Affinity purified

Applications

Tested Applications:

WB, IF, ELISA

Species Specificity:

human

Background Information

BIK, also named as NBK, BIP1 and BP4, is a pro-apoptotic BCL-2 family member and is the founding member of a subfamily of pro-apoptotic proteins known as BH3-alone proteins. It is a transcript that is upregulated during induction of apoptosis.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

