

For Research Use Only

ZNF131 Polyclonal antibody, PBS Only

Catalog Number: 14401-1-PBS



Basic Information

Catalog Number:

14401-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5585

GenBank Accession Number:

BC035875

GeneID (NCBI):

7690

UNIPROT ID:

P52739

Full Name:

zinc finger protein 131

Calculated MW:

71 kDa

Observed MW:

80-85 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IP, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 14401-1-AP (ZNF131 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 14401-1-PBS in a different storage buffer formulation.



IP result of anti-ZNF131 (IP:14401-1-AP, 4ug; Detection:14401-1-AP 1:1000) with HGC-27 cells lysate 1520 ug. This data was developed using the same antibody clone with 14401-1-PBS in a different storage buffer formulation.