

For Research Use Only

beta-actin-Specific Polyclonal antibody



Catalog Number: 14439-1-AP

6 Publications

Basic Information

Catalog Number:

14439-1-AP

Size:

150ul , Concentration: 67 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001101

GeneID (NCBI):

60

UNIPROT ID:

P60709

Full Name:

actin, beta

Calculated MW:

42 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat, zebrafish

Cited Species:

human

Background Information

This antibody is specific to B-actin, but not other actins.

Notable Publications

Author	Pubmed ID	Journal	Application
P Luo	29740907	Neurogastroenterol Motil	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

