

For Research Use Only

TMEM54 Polyclonal antibody

Catalog Number: 14618-1-AP



Basic Information

Catalog Number:

14618-1-AP

Size:

150ul , Concentration: 200 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6234

GenBank Accession Number:

BC001418

GeneID (NCBI):

113452

UNIPROT ID:

Q969K7

Full Name:

transmembrane protein 54

Calculated MW:

24 kDa

Observed MW:

24 kDa, 43 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

TMEM54, also named as BCLP and CAC1, is widely expressed in cancer cell lines.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

