

For Research Use Only

TUBGCP5 Polyclonal antibody, PBS Only

Catalog Number: 14620-1-PBS



Basic Information

Catalog Number:

14620-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6248

GenBank Accession Number:

BC071560

GeneID (NCBI):

114791

UNIPROT ID:

Q96RT8

Full Name:

tubulin, gamma complex associated protein 5

Calculated MW:

118 kDa

Observed MW:

118-120 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 14620-1-AP (TUBGCP5 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 14620-1-PBS in a different storage buffer formulation.



mouse brain tissue were subjected to SDS PAGE followed by western blot with 14620-1-AP (TUBGCP5 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 14620-1-PBS in a different storage buffer formulation.