

For Research Use Only

DNMT3A Polyclonal antibody

Catalog Number: 14623-1-AP **2 Publications**



Basic Information

Catalog Number: 14623-1-AP	GenBank Accession Number: BC043617	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 127 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 1788	
Source: Rabbit	UNIPROT ID: Q9Y6K1	
Isotype: IgG	Full Name: DNA (cytosine-5-)-methyltransferase 3 alpha	
Immunogen Catalog Number: AG6271	Calculated MW: 912 aa, 102 kDa	

Applications

Tested Applications:
ELISA

Cited Applications:
WB

Species Specificity:
human, mouse, rat

Cited Species:
human

Notable Publications

Author	Pubmed ID	Journal	Application
Yanna Ma	35666085	J Chemother	WB
Dongjiao Chen	40975192	Free Radic Biol Med	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data