

For Research Use Only

WIBG Polyclonal antibody

Catalog Number: 14658-1-AP



Basic Information

Catalog Number: 14658-1-AP	GenBank Accession Number: BC006135	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 800 ug/ml by Nanodrop;	GeneID (NCBI): 84305	Recommended Dilutions: IHC 1:50-1:500
Source: Rabbit	UNIPROT ID: Q9BRP8	IF/ICC 1:200-1:800
Isotype: IgG	Full Name: within bgcn homolog (Drosophila)	
Immunogen Catalog Number: AG6276	Calculated MW: 23 kDa	

Applications

Tested Applications: IHC, IF/ICC, ELISA	Positive Controls:
Species Specificity: human, mouse, rat	IHC : mouse colon tissue,
Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0	IF/ICC : U2OS cells,

Background Information

Protein WIBG homolog, also known as PYM1 (Partner of Y14 and mago), is a cytoplasmic RNA-binding protein, that can interact with the ribosome and EJC, may act as an EJC destabilizing factor, preventing EJC accumulation at positions away from exon-exon junctions (non-canonical position).

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded mouse colon tissue slide using 14658-1-AP (WIBG antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunofluorescent analysis of (4% PFA) fixed U2OS cells using WIBG antibody (14658-1-AP) at dilution of 1:400 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L), CL594-phalloidin (red).