

For Research Use Only

LPCAT3 Polyclonal antibody

Catalog Number: 14668-1-AP

3 Publications



Basic Information

Catalog Number:

14668-1-AP

Size:

150ul , Concentration: 180 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6329

GenBank Accession Number:

BC065194

GeneID (NCBI):

10162

UNIPROT ID:

Q6P1A2

Full Name:

lysophosphatidylcholine acyltransferase 3

Calculated MW:

56 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IHC

Species Specificity:

human, mouse, rat

Cited Species:

mouse

Notable Publications

Author	Pubmed ID	Journal	Application
Lijun Zou	39827610	Ecotoxicol Environ Saf	WB
Furui Feng	39260333	Biochem Biophys Res Commun	WB
Ling Kan Chi	38882385	Heliyon	WB,IHC

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

