

For Research Use Only

Muskelin Polyclonal antibody, PBS Only

Catalog Number: 14735-1-PBS

Featured Product



Basic Information

Catalog Number:

14735-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6147

GenBank Accession Number:

BC067825

GeneID (NCBI):

4289

UNIPROT ID:

Q9UL63

Full Name:

muskelin 1, intracellular mediator containing kelch motifs

Calculated MW:

735 aa, 85 kDa

Observed MW:

70-75 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 14735-1-AP (Muskelin antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 14735-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human gliomas using 14735-1-AP (Muskelin antibody) at dilution of 1:50 (under 40x lens). This data was developed using the same antibody clone with 14735-1-PBS in a different storage buffer formulation.



Y79 cells were subjected to SDS PAGE followed by western blot with 14735-1-AP (Muskelin antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 14735-1-PBS in a different storage buffer formulation.