

For Research Use Only

PWP1 Polyclonal antibody

Catalog Number: 14767-1-AP



Basic Information

Catalog Number:

14767-1-AP

Size:

150ul , Concentration: 200 ug/ml by Nanodrop and 133 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6533

GenBank Accession Number:

BC001652

GeneID (NCBI):

11137

UNIPROT ID:

Q13610

Full Name:

PWP1 homolog (S. cerevisiae)

Calculated MW:

56 kDa

Observed MW:

50-56 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : Jurkat cells, human liver tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Jurkat cells were subjected to SDS PAGE followed by western blot with 14767-1-AP (PWP1 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.