

For Research Use Only

LAGE-1 Polyclonal antibody

Catalog Number: 14791-1-AP



Basic Information

Catalog Number:

14791-1-AP

Size:

150ul , Concentration: 340 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6456

GenBank Accession Number:

BC002833

GeneID (NCBI):

30848

UNIPROT ID:

O75638

Full Name:

cancer/testis antigen 2

Calculated MW:

21 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, ELISA

Species Specificity:

human, mouse

Background Information

LAGE-1, also named as ESO2, CT6.2, NY-ESO-2 and CT2, is a potential targets for immunotherapy of some cancer.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data