

For Research Use Only

HMHA1 Polyclonal antibody, PBS Only

Catalog Number: 14832-1-PBS



Basic Information

Catalog Number: 14832-1-PBS	GenBank Accession Number: BC065223	Purification Method: Antigen affinity purification
Size: 100ug , Concentration: 1mg/ml by Nanodrop;	GeneID (NCBI): 23526	
Source: Rabbit	UNIPROT ID: Q92619	
Isotype: IgG	Full Name: histocompatibility (minor) HA-1	
Immunogen Catalog Number: AG6636	Calculated MW: 1136 aa, 125 kDa	
	Observed MW: 140-150 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse, rat

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data