

For Research Use Only

USP10 Polyclonal antibody

Catalog Number: 15025-1-AP



Basic Information

| | | |
|--|--|--|
| Catalog Number: 15025-1-AP | GenBank Accession Number: BC000263 | Purification Method: Antigen affinity purification |
| Size: 150ul , Concentration: 487 µg/ml by Bradford method using BSA as the standard; | GeneID (NCBI): 9100 | |
| Source: Rabbit | UNIPROT ID: Q14694 | |
| Isotype: IgG | Full Name: ubiquitin specific peptidase 10 | |
| Immunogen Catalog Number: AG7106 | Calculated MW: 87 kDa, 92 kDa | |

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse, rat

Background Information

USP10, also named as USP32 and NY-REN-60, belongs to the peptidase C19 family. The deubiquitinating protease USP10 has been identified as a new regulator of p53 in the DNA damage response and tumor development.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data