

For Research Use Only

PTIP Polyclonal antibody

Catalog Number: 15034-1-AP

2 Publications



Basic Information

Catalog Number:

15034-1-AP

Size:

150ul , Concentration: 273 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4567

GenBank Accession Number:

BC033781

GeneID (NCBI):

22976

UNIPROT ID:

Q6ZW49

Full Name:

PAX interacting (with transcription-activation domain) protein 1

Calculated MW:

118 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IP, IF

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse

Notable Publications

Author	Pubmed ID	Journal	Application
Qifan Wang	35354042	Cell Rep	WB,IP
Yan Huang	31220795	Biomaterials	WB,IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

