

For Research Use Only

DYNLRB1 Polyclonal antibody

Catalog Number: 15049-1-AP

3 Publications



Basic Information

Catalog Number:

15049-1-AP

Size:

150ul, Concentration: 200 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7032

GenBank Accession Number:

BC002481

GeneID (NCBI):

83658

UNIPROT ID:

Q9NP97

Full Name:

dynein, light chain, roadblock-type 1

Calculated MW:

11 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IF, ELISA

Cited Applications:

IF

Species Specificity:

human, mouse, rat

Cited Species:

human, rat

Notable Publications

Author	Pubmed ID	Journal	Application
Syeda Ridita Sharif	27646688	Mol Cells	IF
Md Ariful Islam	26467288	Mol Cells	IF
Md Ariful Islam	26272270	Exp Mol Med	IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

