

For Research Use Only

G Alpha 13 Polyclonal antibody

Catalog Number: 15142-1-AP



Basic Information

Catalog Number:

15142-1-AP

Size:

150ul , Concentration: 1260 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7316

GenBank Accession Number:

BC036756

GeneID (NCBI):

10672

UNIPROT ID:

Q14344

Full Name:

guanine nucleotide binding protein (G protein), alpha 13

Calculated MW:

44 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human testis tissue were subjected to SDS PAGE followed by western blot with 15142-1-AP (G alpha 13 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.