

For Research Use Only

SMPD2 Polyclonal antibody

Catalog Number: 15290-1-AP

1 Publications



Basic Information

Catalog Number:

15290-1-AP

Size:

150ul , Concentration: 107 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7482

GenBank Accession Number:

BC000038

GeneID (NCBI):

6610

UNIPROT ID:

O60906

Full Name:

sphingomyelin phosphodiesterase 2, neutral membrane (neutral sphingomyelinase)

Calculated MW:

48 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

IF

Species Specificity:

human, mouse, rat

Cited Species:

mouse

Notable Publications

Author	Pubmed ID	Journal	Application
Wei Fan	34042046	Elife	IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

