

For Research Use Only

# TRIAP1 Polyclonal antibody

Catalog Number: 15351-1-AP

4 Publications



## Basic Information

Catalog Number:

15351-1-AP

Size:

150ul, Concentration: 1000 ug/ml by Nanodrop and 380 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7597

GenBank Accession Number:

BC002638

GeneID (NCBI):

51499

UNIPROT ID:

O43715

Full Name:

TP53 regulated inhibitor of apoptosis

1

Calculated MW:

9 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IHC 1:50-1:500

## Applications

Tested Applications:

IHC, ELISA

Cited Applications:

WB, IHC

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse

**Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (\*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0**

Positive Controls:

IHC : mouse pancreas tissue, human breast cancer tissue

## Notable Publications

Author	Pubmed ID	Journal	Application
Julia Ketteler	30871022	J Clin Med	WB, IHC
Guohua Liu	39135128	Clin Respir J	WB
Zhi-Ying Shao	38769122	Cell Death Dis	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunohistochemical analysis of paraffin-embedded mouse pancreas tissue slide using 15351-1-AP (TRIAP1 antibody) at dilution of 1:200 (under 10x lens)..



Immunohistochemical analysis of paraffin-embedded mouse pancreas tissue slide using 15351-1-AP (TRIAP1 antibody) at dilution of 1:200 (under 40x lens)..