

For Research Use Only

AAK1 Polyclonal antibody

Catalog Number: 15353-1-AP



Basic Information

Catalog Number:

15353-1-AP

Size:

150ul

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7496

GenBank Accession Number:

BC002695

GeneID (NCBI):

22848

UNIPROT ID:

Q2M2I8

Full Name:

AP2 associated kinase 1

Calculated MW:

104 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

AAK1, also named as KIAA1048, belongs to the protein kinase superfamily and Ser/Thr protein kinase family. It phosphorylates the AP2M1/mu2 subunit of the adaptor protein complex 2 (AP-2). AAK1 may play a role in regulating aspects of clathrin-mediated endocytosis (PMID:14617351).

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

