

For Research Use Only

CDK2AP2 Polyclonal antibody

Catalog Number: 15423-1-AP



Basic Information

Catalog Number:

15423-1-AP

Size:

150ul , Concentration: 700 ug/ml by Nanodrop and 333 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7730

GenBank Accession Number:

BC002850

GeneID (NCBI):

10263

UNIPROT ID:

O75956

Full Name:

cyclin-dependent kinase 2 associated protein 2

Calculated MW:

13 kDa

Observed MW:

16 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HeLa cells,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 15423-1-AP (CDK2AP2 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.