For Research Use Only

NOL10 Polyclonal antibody

Catalog Number: 15433-1-AP



Basic Information

Catalog Number:

GenBank Accession Number:

Purification Method:

15433-1-AP

Size:

BC005125 GeneID (NCBI): Antigen affinity purification

150ul , Concentration: 173 µg/ml by 79954 Bradford method using BSA as the standard;

UNIPROT ID: Q9BSC4 Full Name:

Source: Rabbit Isotype:

nucleolar protein 10 Calculated MW:

IgG Immunogen Catalog Number:

AG7692

Applications

Tested Applications:

ELISA

Species Specificity: human, mouse, rat

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1%BSA

Selected Validation Data