

For Research Use Only

# TRADD Polyclonal antibody, PBS Only

Catalog Number: 15468-1-PBS



## Basic Information

Catalog Number:

15468-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7768

GenBank Accession Number:

BC004491

GeneID (NCBI):

8717

UNIPROT ID:

Q15628

Full Name:

TNFRSF1A-associated via death domain

Calculated MW:

34 kDa

Observed MW:

34 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, IP, Indirect ELISA

Species Specificity:

human, mouse

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 15468-1-AP (TRADD antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 15468-1-PBS in a different storage buffer formulation.



IP result of anti-TRADD (IP:15468-1-AP, 4ug; Detection:15468-1-AP 1:500) with Jurkat cells lysate 4000 ug. This data was developed using the same antibody clone with 15468-1-PBS in a different storage buffer formulation.