

For Research Use Only

# FOXP3 Polyclonal antibody

Catalog Number: 15496-1-AP



## Basic Information

Catalog Number:	GenBank Accession Number:	Purification Method:
15496-1-AP	NM_014009	Antigen affinity purification
Size:	GeneID (NCBI):	
150ul	50943	
Source:	Full Name:	
Rabbit	forkhead box P3	
Isotype:	Calculated MW:	
IgG	47 kDa	

## Applications

Tested Applications:  
ELISA

Species Specificity:  
human, mouse, rat

## Background Information

FOXP3, also named as IPEX, JM2 and Scurfin, is a transcription factor. FOXP3 plays a critical role in the control of immune response. Defects in FOXP3 are the cause of immunodeficiency polyendocrinopathy, enteropathy, X-linked syndrome (IPEX) which also known as X-linked autoimmunity-immunodeficiency syndrome. This antibody is a rabbit polyclonal antibody raised against a peptide mapping within human FOXP3.

## Storage

Storage:  
Store at -20°C. Stable for one year after shipment.

Storage Buffer:  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

