

For Research Use Only

DISC1 Polyclonal antibody

Catalog Number: 15501-1-AP



Basic Information

Catalog Number:

15501-1-AP

Size:

150ul , Concentration: 167 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001012959

GeneID (NCBI):

27185

UNIPROT ID:

Q9NRI5

Full Name:

disrupted in schizophrenia 1

Calculated MW:

94 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

Disrupted-in-schizophrenia 1 (DISC1) is a genetic susceptibility factor for schizophrenia and related severe psychiatric conditions, with several translation isoforms. DISC1 isoforms have distinct developmental patterns of expression. This antibody binds both short and long isoforms of DISC1.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data