

For Research Use Only

SSBP4 Polyclonal antibody

Catalog Number: 15595-1-AP



Basic Information

Catalog Number:

15595-1-AP

Size:

150ul , Concentration: 400 ug/ml by Nanodrop and 347 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7952

GenBank Accession Number:

BC000274

GeneID (NCBI):

170463

UNIPROT ID:

Q9BWG4

Full Name:

single stranded DNA binding protein 4

Calculated MW:

39 kDa

Observed MW:

39-42 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : U-937 cells,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



U-937 cells were subjected to SDS PAGE followed by western blot with 15595-1-AP (SSBP4 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.